


The Promise of PFA Delivered

Demonstrated safety and clinical excellence in one complete system



PULSAR-IDE

Trial summary

The global, multicenter PULSAR IDE Study generated data supporting the safety and effectiveness of the Globe® Pulsed Field System in treating patients with paroxysmal atrial fibrillation.



12 Sites



19 Operators



4 Countries

United States, Canada,
Germany, Czechia



183 Patients

Trial design

- Prospective, multicenter, single-arm Investigational Device Exemption trial
- 183 total patients with symptomatic paroxysmal atrial fibrillation enrolled
- 164 of 183 patients included in the safety and effectiveness analysis, with 19 as lead-in cohort
- All patients underwent pulmonary vein isolation only

12-month follow-up

- Weekly and symptomatic transtelephonic monitoring
- 12-lead ECG at 3, 6, and 12 months
- 24-hour Holter monitoring at 6 and 12 months

At a glance

0.6%

Primary Safety
Event Rate

0%

Device-Related Primary
Safety Event Rate

77.8%

12-Month Primary
Effectiveness, Paroxysmal Cohort

Safety

No device-related primary adverse safety events reported

- Zero reports of cardiac tamponade, atrio-esophageal fistula, or phrenic nerve injury
- One non-device-related incident of a hemorrhagic stroke following ablation in a patient with severe hypertension

0/164 Death
0/164 Cardiac tamponade/perforation
0/164 Heart block
0/164 Major vascular access complication
0/164 Major bleeding
0/164 Myocardial infarction (MI)
0/164 Pericarditis requiring intervention
0/164 Pulmonary edema
0/164 Transient ischemic attack (TIA)
0/164 Vagal nerve injury resulting in gastroparesis
0/164 Stoke/cerebrovascular accident
0/164 Atrioesophageal fistula*
0/164 Diaphragmatic paralysis/phrenic nerve injury†
0/164 Severe PV stenosis‡
0/164 Any Globe PF System-related or PF procedure-related cardiovascular and/or pulmonary adverse event within 7 days that prolongs or requires hospitalization for more than 48 hours§

* within 2 months of ablation procedure

† unresolved at 6 months

‡ within 6 months of ablation procedure

§ excluding recurrent AF/AFL/AT

¶ n=10


Effectiveness

77.8% freedom from atrial arrhythmias at 1 year

The primary effectiveness endpoint was defined as freedom from AF/AFL/AT at 12 months, acute procedure failure, direct current cardioversion, antiarrhythmic drug failure, redo-ablation due to recurrent AF/AFL/AT, and surgical treatment for AF/AFL/AT.


Durable isolation

- Only 6% of patients required repeat ablation¶
- 95% (35/37) of assessed veins remained durably isolated upon repeat procedures¶




100% Acute Procedural Success*

* Defined as acute isolation of all accessible targeted PV, confirmed by entrance and/or exit block with the Globe Catheter.




Follow-Up and Compliance

- Completed 12-month follow-up compliance of **99%**
- Weekly and symptomatic transtelephonic monitoring (TTM) with a compliance of **71%**
- 12-lead ECG at 3, 6, and 12 months with a compliance of **95%**
- 24-hour Holter monitoring at 6 and 12 months with a compliance of **91%**




Live local array

provides localized beat-by-beat data on region of interest




CONTACT™ Scan

confirms electrode-tissue contact for therapy planning




Tailored electrode selection,
ensures only the intended tissue is targeted



Efficiency

The Globe Pulsed Field System provides true single-shot ablation, averaging one application per vein.




Average of 1.18 pulsed field applications per vein to achieve successful acute isolation of the targeted PV



See the anatomy. See the target.
See the contact.

With fully integrated HD mapping, the Globe System creates fast, real-time visuals of atrial anatomy and activity.



Characteristics	PVI only (n = 164)
Procedure time, min	95.8 ± 20.7
Left atrial dwell time, min	59.9 ± 11.6
Transpired ablation time, min	25.5 ± 9.9
PFA delivery time, s	49.2 ± 10.2
Fluoroscopy time, min	9.2 ± 7.5

Values are presented as mean ± SD.
PFA = pulsed field ablation; PVI = pulmonary vein isolation.

The Globe® PF System is investigational and not approved for sale in any geography.
CAUTION - Investigational device. Limited by Federal (or United States) law to investigational use.



Kardium Inc.

155 – 8518 Glenlyon Parkway | Burnaby, BC V5J 0B6 | CANADA

kardium.com

©2025 Kardium Inc. All rights reserved. All trademarks are the property of their respective owners.